How to use

BD Veritor™ tests

Detect COVID-19



COLLECTION





Roll **5 times** in **each nostril**; right and left nostrils

2

PREPARATION



Remove cap from tube and insert nasal swab



Plunge swab up and down 15 seconds




Pinch tube slightly to extract liquid while removing swab



Firmly place cap on tube




TESTING



Add **3 drops** to the test strip sample well

Wait 15 min.
Here is an
example of
a positive test and
a negative test.



If you are positive, please call the Info-Covid line: 1-888-662-7482

